

# RESIST- 4 O.K.N.V.



www.corisbio.com  
IFU- 58R8/TB/01

Manufacturer:

Coris BioConcept  
Science Park CREALYS  
Rue Jean Sonet 4A  
B – 5032 GEMBLoux  
BELGIUM  
Tel.: +32(0)81.719.917  
Fax: +32(0)81.719.919  
info@corisbio.com

Produced in BELGIUM

## Rapid diagnostic test for the detection of OXA-48, KPC, NDM and VIM carbapenemases on bacterial colony

FOR *IN VITRO* DIAGNOSTIC USE  
FOR PROFESSIONAL USE ONLY

EN

References: K-15R8, 2x20 cassettes, buffer, 20 tubes and droppers

(EN) For Instructions For Use in your language : (FR) Pour obtenir les notices dans la langue de votre choix : (ES) Para las instrucciones de uso en su idioma : (PT) Para Instruções de Uso na sua língua : (IT) Per le Istruzioni di Uso nella sua lingua : (DE) Für Gebrauchsanleitungen in Ihrer Sprache : (NL) Voor Gebruiksaanwijzing in uw eigen taal :	
	www.e-labeling.eu/cor58R8
	(EU) +800 135 79 135 (non-EU) +31 20 794 7071 (CA) +1 855 805 8539 (AR, CO, UY, AU, NZ) +800 135 79 135

### I. INTRODUCTION

Carbapenemase-producing organisms (CPO) and more particularly carbapenem-resistant Enterobacteriaceae (CRE) represent a major public health concern worldwide due to their broad spectrum of resistance to antibiotics including, besides carbapenems, most classes of antimicrobial agents, and thus leaving very few options for the management of infected patients. The rapid spread of CPOs or of the genes encoding these resistances has led to nosocomial outbreaks and endemic situations in several countries in Europe and elsewhere worldwide.

Development of new rapid diagnostic tests to track antimicrobial resistance patterns is considered as one of the priority core action by international experts and health authorities. NDM and KPC represent two of the most increasing and prevalent carbapenemases in many countries. On the other hand, class D OXA-48 type carbapenemases represent the most challenging resistance mechanism to detect for clinical laboratories. VIM is not only present in Enterobacteriaceae but is also highly prevalent in non-fermenter bacteria. Rapid identification of those carbapenemases is of utmost importance to improve both patient therapy and control of the spread of such antibiotic resistance in hospitals.

Confirmatory phenotypic tests using combination disks with specific inhibitors already exist for detection of selected types of carbapenemases including class A (KPC) and class B (VIM, IMP, NDM) carbapenemases; however, these tests are time-consuming and require an extra additional day following antimicrobial susceptibility testing results. Moreover, phenotypic colorimetric assays are in some instances not sensitive enough for the detection of low-activity carbapenemases such as OXA-48. Several molecular assays based on different formats also allow detection of carbapenemases. These tests are expensive, time-consuming and can only be performed in dedicated environment and by skilled personnel, hence limiting their generalized use.

### II. PRINCIPLE OF THE TESTS

These tests are ready to use and are based on a membrane technology with colloidal gold nanoparticles. Each pouch contains 2 lateral-flow cassettes for the identification of (i) KPC and OXA-48 and (ii) NDM and VIM. These two devices are aimed at the detection of KPC, OXA-48, NDM and VIM carbapenemases on a single colony of bacterial isolates growing on agar plate resuspended in the provided buffer.

**Identification of KPC and OXA-48.** A nitrocellulose membrane is sensitized with:

- (1) a monoclonal antibody directed against the KPC carbapenemase (bottom K line)
- (2) a monoclonal antibody directed against the OXA-48 carbapenemase (middle O line)
- (3) a control capture reagent (upper C line).

Three different colloidal gold nanoparticles conjugates are dried on a membrane: a conjugate directed against the KPC carbapenemase, a conjugate directed against the OXA-48 carbapenemase, and a control conjugate.

**Identification of NDM and VIM.** A nitrocellulose membrane is sensitized with:

- (1) a monoclonal antibody directed against the NDM carbapenemase (bottom N line),
- (2) a monoclonal antibody directed against the VIM carbapenemase (middle V line),
- (3) a control capture reagent (upper C line).

Three different colloidal gold nanoparticles conjugates are dried on a membrane: a conjugate directed against the NDM carbapenemase, a conjugate directed against the VIM carbapenemase, and a control conjugate.

When the provided buffer containing the resuspended bacteria comes into contact with the strip, the solubilised conjugates migrate with the sample by passive diffusion and conjugates and sample material come into contact with the immobilized respective antibodies that are adsorbed onto the nitrocellulose strip. If the sample contains a KPC, OXA-48, NDM or VIM carbapenemase, the respective complexes made of the conjugates and either KPC, or OXA-48, or NDM or VIM will remain bound to their respective specific lines (KPC, K line; OXA-48, O line; NDM, N line, VIM, V line). The migration continues by passive diffusion and both conjugates and sample material come into contact with the (upper) line control reagent that binds a control conjugate (line C), thereby producing a red line. The result is visible within 15 minutes in the form of red lines on the strip.

### III. REAGENTS AND MATERIALS

#### 1. RESIST-4 O.K.N.V. (2x20 cassettes)

20 sealed pouches containing two lateral-flow cassettes and one desiccant.

Each device contains one sensitized strip.

#### 2. LY-A buffer vial (15 mL)

Saline solution buffered to pH 7.5 containing TRIS, Na<sub>3</sub> (<0,1%) and a detergent.

#### 3. Instruction for use (1)

#### 4. Semi-rigid disposable collection tubes with droppers (20)

### IV. SPECIAL PRECAUTIONS

All operations linked to the use of the test must be performed in accordance with Good Laboratory Practices (GLP).

All reagents are for *in vitro* diagnostic use only.

Pouch must be opened with care:

Avoid touching nitrocellulose with your fingers.

Wear gloves when handling samples.

Never use reagents from another kit.

Green lines indicate immunoreagents adsorption sites. Green colour disappears during the test.

Reagents' quality cannot be guaranteed beyond their shelf-life dates or if reagents are not stored under required conditions as indicated in the insert.

### V. WASTE DISPOSAL

Dispose of gloves, swabs, test tubes and used devices in accordance with GLP.

Each user is responsible for the management of any waste produced, and must ensure that it is disposed of in accordance with the applicable legislation.

### VI. STORAGE

An unopened pouch may be kept at between 4 and 30°C and used until the shelf-life date indicated on the packaging. Once the pouch is opened, run the test immediately.

Avoid freezing devices and buffer.

### VII. SPECIMEN HANDLING AND COLLECTION

Specimens to be tested should be obtained and handled by standard microbiological methods.

Make sure that the specimens are not treated with solutions containing formaldehyde or its derivatives.

### VIII. PROCEDURE

PREPARATIONS OF THE TEST:


Allow kit components, in unopened packaging, and specimens (in case the plate containing colony to be tested was kept at 4°C) to reach room temperature (15-30°C) before performing a test.

Open the pouch and remove the device. Once opened, run the test immediately. Indicate the patient's name or specimen number on the device (one device per sample).

SPECIMEN PREPARATION PROCEDURE:

Performance claims with regard to samples types other than bacterial colonies have not been established. We recommend the use of fresh bacterial colonies for optimal test performance.

1. Prepare one semi-rigid tube and add 12 drops of LY-A buffer in the tube.
2. Harvest bacteria by touching one colony with a disposable bacteriological loop and dip the loop till the bottom of the semi-rigid tube containing the buffer.
3. Stir thoroughly to homogenize the solution.
4. Insert tightly the dropper on the semi-rigid tube.
5. Invert the test tube and add slowly 3 drops of diluted sample into the sample well of each of the two cassettes labeled (i) KPC and OXA-48 and (ii) NDM and VIM.
6. Allow to react for 15 min max and read the result.



Positive results may be reported as soon as the test and control lines become visible.

Do not take the appearance of new lines into account after the reaction time is passed.

The result must be read on still wet strip.

### IX. INTERPRETING RESULTS

The results are to be interpreted as follows for each of the two cassettes:

**Negative test result:** a reddish-purple line appears across the central reading window at the Control line (C) position. No other band is present.


**Positive test result:** in addition to a reddish-purple band at the Control line (C), a visible reddish-purple band appears at one of the Test lines position (OXA-48 or KPC) on cassette labelled (i) KPC and OXA-48, or at one of the Test lines position (VIM or NDM) on cassette labelled (ii) NDM and VIM. Intensity of the test line may vary according to the quantity of antigens as well as of the variant type present in the sample. Any reddish-purple Test line (OXA-48, KPC, NDM and VIM), even weak, should be considered as a positive result.

If a positive test line appears beside of the O mark, the sample contains OXA-48 or OXA-48-like variant, beside of K mark, the sample contains KPC, beside of N mark, the sample contains NDM and beside of V mark, VIM is present in the sample. Combinations of positive test lines can occur. In this case the sample contains the combination of several carbapenemases.


**Invalid test result:** The absence of a Control line indicates a failure in the test procedure. Repeat invalid tests with a new test device.

Note: during the drying process, a very faint shadow may appear at the Test line positions (O, K, N, V). It should not be regarded as a positive result.

Cassette 1 : OXA-48 & KPC



Cassette 2 : VIM & NDM



**X. PERFORMANCE**

**A. Detection Limit**

The detection limit determined with purified recombinant proteins of OXA-48, KPC, NDM and VIM have been evaluated at 0.125 ng/ml, 0.625 ng/ml, 0.25 ng/ml and 0.23 ng/ml, respectively.

**B. Prospective study (based on RESIST-3 O.K.N. kit)**

The OXA-48 and KPC cassette test was validated by comparison with reference molecular method (validated multiplex PCR including sequencing) in the National Reference Laboratory for Multidrug-Resistant Gram Negative Bacilli (Belgium) in a prospective study performed on 173 non duplicated, consecutive suspected CPE clinical isolates referred from July to September 2016.

OXA-48 test	Molecular method	Positive	Negative	Total
	<b>Positive</b>	69	0	69
	<b>Negative</b>	0	104	104
	<b>Total</b>	69	104	173

95 % Confidence Interval <sup>1</sup>

Sensitivity:	100 %	(95.7 to 100 %)
Specificity:	100 %	(97.2 to 100 %)
Positive Predictive value:	100 %	(95.7 to 100 %)
Negative predictive value:	100 %	(97.2 to 100 %)
Agreement:	100 %	(173/173)

KPC test	Molecular method	Positive	Negative	Total
	<b>Positive</b>	9	0	9
	<b>Negative</b>	0	164	164
	<b>Total</b>	9	164	173

95 % Confidence Interval <sup>1</sup>

Sensitivity:	100 %	(68.4 to 100 %)
Specificity:	100 %	(98.2 to 100 %)
Positive Predictive value:	100 %	(68.4 to 100 %)
Negative predictive value:	100 %	(98.2 to 100 %)
Agreement:	100 %	(173/173)

**C. Validation on collection of reference strains**

The VIM and NDM cassette test was validated by comparison with reference molecular method in the National Reference Laboratory for Multidrug-Resistant Gram Negative Bacilli (Belgium) in a retrospective study.

NDM test	Molecular method	Positive	Negative	Total
	<b>Positive</b>	24	0	24
	<b>Negative</b>	0	95	95
	<b>Total</b>	24	95	119

95 % Confidence Interval <sup>1</sup>

Sensitivity:	100 %	(82.8 to 100 %)
Specificity:	100 %	(95.2 to 100 %)
Positive Predictive value:	100 %	(82.8 to 100 %)
Negative predictive value:	100 %	(95.2 to 100 %)
Agreement:	100 %	(119/119)

VIM test	Molecular method	Positive	Negative	Total
	<b>Positive</b>	38	0	38
	<b>Negative</b>	1*	80	81
	<b>Total</b>	39	80	119

\*: the false-negative result is a *P. aeruginosa* colony harboring VIM-5 and NDM-1 genes. This colony was detected as NDM-positive but VIM-negative. The production of VIM-5 was not assessed.

Sensitivity:	97.4 %	95 % Confidence Interval <sup>1</sup>
Specificity:	100 %	(84.9 to 99.9 %)
Positive Predictive value:	100 %	(94.3 to 100 %)
Negative predictive value:	98.8 %	(88.6 to 100 %)
Agreement:	99.2 %	(92.4 to 99.9 %)
		(118/119)

**D. Repeatability and reproducibility**

To check intra-batch accuracy (repeatability), the same positive samples and a buffer solution were processed 15 times on kits of the same production batch in the same experimental conditions. All observed results were confirmed as expected. To check inter-batch accuracy (reproducibility), some samples (positive and buffer) were processed on kits from three different production batches. All results were confirmed as expected.

**XI. LIMITS OF THE KIT**

The test is qualitative and cannot predict the quantity of antigens present in the sample. Clinical presentation and other test results must be taken into consideration to establish diagnosis.

A positive test does not rule out the possibility that other antibiotic resistance mechanisms may be present.

**XII. TECHNICAL PROBLEMS / COMPLAINTS**

If you encounter a technical problem:

1. If possible, keep the sample in the appropriate storage condition during the complaint management
2. Contact Coris BioConcept ([client.care@corisbio.com](mailto:client.care@corisbio.com))

**XIII. BIBLIOGRAPHIC REFERENCES**

A. Y. Glupczynski, A. Jousset, S. Evrard, R. Bonnin, T. Huang, L. Dortet, P. Bogaerts and T. Naas. Prospective evaluation of the OKN K-SeT assay, a new multiplex immunochromatographic test for the rapid detection of OXA-48-like, KPC and NDM carbapenemases. J Antimicrob Chemother. 2017 Mar (7):1955-1960

B. DW. Wareham and MH. Abdul Momin Rapid Detection of Carbapenemases in Enterobacteriaceae: 1 Evaluation of the RESIST-3 O.K.N (OXA-48, KPC, NDM) Multiplexed Lateral Flow Assay. J Clin Microbiol. 2017 Feb (4):1223-1225

C. F Pasteran, L Denorme, I Ote, S Gomez, D De Belder, Y Glupczynski, P Bogaerts, B Ghiglione, P Power,

D. P Mertens, A Corso. Rapid identification of OXA-48 and OXA-163 subfamily in carbapenem resistant gram-negative bacilli with a novel immunochromatographic lateral flow assay. J Clin Microbiol. 2016 Aug 17. pii: JCM.01175-16

E. D. Meunier, A. Vickers, R. Pike, R.L. Hill, N. Woodford and K.L. Hopkins. Evaluation of the K-SeT R.E.S.I.S.T. immunochromatographic assay for the rapid detection of KPC and OXA-48-like carbapenemases. J Antimicrob Chemother. 2016 Aug; 71 (8):2357-9

F. Wareham DW, Shah R, Betts JW, Phee LM, Abdul Momin MH. Evaluation of an Immunochromatographic Lateral Flow Assay (OXA-48 K-SeT) for the Rapid Detection of OXA-48-like Carbapenemases in Enterobacteriaceae. J Clin Microbiol. 2016 Feb;54 (2):471-3

G. Fernández J, Fleites A, Rodicio MR, Vazquez F. Evaluation of OXA-48 K-Se T: an immunochromatographic assay for rapid detection of OXA-48-producing Enterobacteriaceae. Diagn Microbiol Infect Dis. 2016 May;85 (1):12-5

I. Dortet L, Jousset A, Sainte-Rose V, Cuzon G, Naas T. Prospective evaluation of the OXA-48 K-SeT assay, an immunochromatographic test for the rapid detection of OXA-48-type carbapenemases. J Antimicrob Chemother. 2016 Jul;71 (7):1834-40

J. Glupczynski Y, Evrard S, Ote I, Mertens P, Huang TD, Leclipteux T, Bogaerts P. Evaluation of two new commercial immunochromatographic assays for the rapid detection of OXA-48 and KPC carbapenemases from cultured bacteria. J Antimicrob Chemother. 2016 May;71(5):1217-22

K. B.M. Willey, X. Trimi, R. Ioboni, D.A. Boyd, G. Ricci, D.N. Grohn, D. Terenzi, A. Mazzulli, L. Mataseje, M. Mulvey, P. Lo, T. Mazzulli, S.M. Poutanen. The Coris BioConcept OXA48 K-SeT Immuno-Chromatographic Assay Detects OXA48-type Carbapenemases with High Sensitivity and Specificity. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 9-12, 2016

L. A. Abulaila, F. Erdem, Z. Aktas and O. Oncul Evaluation Comparison of a novel OXA-48 K-Set test and blue-carba test in detection of carbapenemase-producing enterobacteriaceae with using PCR as reference method. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016

M. O. Karatuna, M. Kaya, I. Akyar Evaluation of the performance of OXA-48 K-SeT immunochromatographic test for rapid identification of OXA-48 carbapenemase-producing Enterobacteriaceae. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016

N. C.S. Nodari, A. Barth, C. Magagnin, A. Zavascki, A. Gales, C. Carvalhoes OXA-370 is rapidly detected from different culture media using OXA-48 K-SeT® immunochromatography. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016

P. Bogaerts, S. Evrard, G. Cuzon, TD. Huang, T. Naas and Y. Glupczynski Specificity of the OXA-48 immunochromatographic K-SeT for the detection of OXA-48 like in Shewanella spp. 26th European Congress of Clinical Microbiology and Infectious Diseases, Infectious Amsterdam April 09 – 12, 2016

Q. A. Sarría, R. Gomez-Gil, G. Ruiz, MP. Romero, J. García- Rodríguez Preliminary study of the OXA-48 card test method for the direct detection of OXA-48 carbapenemase in blood and plates culture. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016

R. L. Dortet, A. Jousset, V. Sainte-Rose, G. Cuzon, and T. Naas Prospective evaluation of the OXA-48 K-SeT for the detection of OXA-48-type carbapenemase producing Enterobacteriaceae. 26th European Congress of Clinical Microbiology and Infectious Diseases, Amsterdam April 09 – 12, 2016

Last update: JANUARY 2018

	Catalogue number		Manufacturer
	In vitro diagnostic medical device		Temperature limits
	Contains sufficient for <n> tests		Batch code
	Consult instructions for use		Do not reuse
	Keep dry		Use by
	Diluent specimen		Contains Sodium azide

<sup>1</sup> Newcombe, Robert G. "Two-Sided Confidence Intervals for the Single Proportion: Comparison of Seven Methods," *Statistics in Medicine*, 17, 857-872 (1998).