



Positive Controls

code	name	content	format	dilution	No. tests
C-1099	PYLORI POSITIVE CONTROL	inactivated lysate	liquid (1ml)	dilution 1/2	4 tests
C-1091	E. COLI O157 POSITIVE CONTROL	inactivated lysate	liquid (1ml)	dilution 1/2	4 tests
P-1022	STREP-A POSITIVE CONTROL	inactivated lysate	liquid (0.4ml)	dilution 1/5	4 tests
C-1082	ADENOVIRUS POSITIVE CONTROL	heat-inactivated - type 40	Freeze-dried	dilution 1/2	4 tests
X-1081	ROTAVIRUS POSITIVE CULTURE	culture not inactivated	Freeze-dried	dilution 1/4	8 tests
C-1090	INFLUENZA A POSITIVE CONTROL	recombinant protein	liquid (1ml)	dilution 1/2	4 tests
C-1086	RSV POSITIVE CONTROL	inactivated lysate	liquid (1ml)	dilution 1/2	4 tests
TBC	GIARDIA POSITIVE CONTROL	<i>coming soon</i>			
C1085	CRYPTOSPORIUM POSITIVE CONTROL	inactivated oocysts	Freeze-dried	dilution 1/4	16 tests

Please note

Coris positive controls have been selected for use with Coris kits, however many of these reagents are suitable for use with other manufacturers kits.

If you are unsure about the compatibility between Coris controls and your current assays, please use the contact form to let us know which kits you are using and we can advise you on compatibility issues.